

The Waived Testing Requirements

Employee Name: _____ Date: _____ Score: _____

1. Glucose test strips are affected by humidity; the desiccant in lid prevents strips from becoming defective due to exposure of humidity. It is unacceptable to leave vials uncapped or to leave unused test strips out of their protective, properly capped vial; improperly stored test strips may result in falsely elevated glucose results.

An indication that test strips have been damaged (i.e. by humidity) would be reflected in the control results: **Level 1** would be out of range while **Level 2** would be within range.

- a. True
- b. False

2. Use of the Accucheck Inform is limited to screening and monitoring purposes; therefore, a diagnosis of diabetes cannot be made just based on an elevated test result.

- a. True
- b. False

3. Glucose control solutions are only good for _____ after opening. The new expiration date for a brand new control set opened on 1/1/07 whose manufacturer expiration date is 3/28/07 is _____, and should be written on both vials to prevent their use beyond expiration.

- a. 6 months; 7/1/06
- b. 3 months; 4/1/07
- c. 3 months; 3/28/07
- d. 9 months; 10/1/07

4. Prior to performing a finger stick, the patient must wash their hands if possible. Then the operator must wipe the fingertip with alcohol and allow their hair to dry. Failure to thoroughly clean the fingertip and eliminate any sugary food remnants (for example, after the patient has eaten breakfast) may lead to these substances being measured as part of the patient's blood sugar or a falsely elevated glucose test.

- a. True
- b. False

5. Maintenance of the Accucheck Inform requires wiping of meter surfaces using a cloth dampened with _____. Whichever agent is used, one should never spray meter directly, or submerge it in solution. Also, the meter must never be exposed to high temperatures (i.e. isolettes) for longer than _____.

- a. 10% bleach or 70% isopropyl alcohol BUT NOT Cavicide because it damages the meter's screen; 10 seconds
- b. Cavicide or 10% bleach BUT NOT 70% isopropyl alcohol because it damages the meter's screen; 10 seconds
- c. Cavicide or 70% isopropyl alcohol BUT NOT 10% bleach because it damages the meter's screen; 10 seconds
- d. I give up...way too confusing, we use cavicide anyways!; 10 seconds

6. Identify the properly dosed test strip by clicking on test strips 1, 2, 3 and 4 below. For improperly dosed test strips, reapplication of sample is acceptable as long as it is done with _____ seconds.
Note: To view test strip image must click on the individual blue link below.
- Test strip 1; 15 seconds
 - Test strip 2; 15 seconds
 - Test strip 3; 20 seconds
 - Test strip 4; 15 minutes
7. A critical blood glucose result for an adult is a value that is _____. Critical results are documented on the Accuchek Inform using the _____ comment code.
- Less than 50 mg/dl or greater than 450 mg/dl; notified provider
 - Less than 10 mg/dl or greater than 600 mg/dl; notified provider
 - Less than 30 mg/dl or greater than 300 mg/dl; lab study
 - Less than 40 mg/dl or greater than 110 mg/dl; operator error
8. As of 2005, Joint Commission requires documentation of blood glucose results in the medical record include:
- The control lot number in use
 - The strip lot number in use
 - The glucometer serial number
 - The glucose reference (normal) range
9. Interfering substances are those substances which may lead to lower or higher BG results (depending on substance type). Manufacturers are obligated to disclose those substances.
As an Accuchek Inform user, you should be aware of the following interfering substances and how each one impacts the BG result.
- Maltose, Xylose, Galactose (elevate result)
 - Estrogens (decrease result)
 - Bilirubin levels > 20 mg/dl (elevate result)
 - a & c
10. Controls must be performed when:
- A new lot of strips is opened, or improper storage is suspected (i.e. refrigerated strips, or vial left open)
 - Once per shift
 - The meter has been dropped
 - The patient's clinical symptom's do not really correlate with the BG result
 - All of the above

NURSING FIELD VISIT COMPETENCY EVALUATION

SKILLED NURSE: _____ POSITION: _____

PATIENT NAME: _____ DATE: _____

AREAS OF EVALUTION

HOME GLUCOSE MONITORING	YES	N/A	NO
Correctly verbalizes purpose of test?	X		
Demonstrates correct procedure for preparing and collecting the sample?	X		
Correctly preserves/handles specimen?	X		
Demonstrates proper use/performance of equipment?	X		
Performs/interprets tests correctly?	X		
Correctly calibrates the instruments?	X		
Implements organization quality control procedures?	X		
Comments: _____ _____ _____			

SUPERVISOR: _____ DATE: _____